



900 Series

COLPOSCOPES

Seiler Precision Microscopes is proud to offer the 900 Series Colposcopes. These colposcopes have state-of-the-art optics, exceptionally bright illumination, and a power supply with (2) ports for the fiber optic cable.

Seiler is a name you can trust for all of your microscopic needs. We have been doing this since 1945. Our well-trained team of microscope consultants can assist you in determining which microscope is best for you. Our service department offers a full-range of repairs and is also available to answer any technical questions that you might have. Take advantage of our great service and very competitive prices and you'll see why Seiler is a name you can trust.



3433 Tree Court Industrial Blvd. • St. Louis, MO 63122
(314) 968-2282 • Toll-Free: (800) 489-2282 • Fax: (314) 968-3601
E-mail: micro@seilerinst.com • www.seilermicro.com



935



PRODUCT #	DESCRIPTION	PRICE
30935-3	Model 935 3-Step Colposcope	\$4,795.00

Includes the following:

- Smooth three-step turret magnification changer.
- 300mm working distance objective lens offering 3x, 7x, and 17x magnifications. Optional magnifications/objective lens combinations are available.
- Widefield 12.5x eyepieces with diopter locks and fold-down rubber eyecups for eyeglass wearers.
- Straight binocular head f=160.
- Fan cooled, **dual-port variable fiber optic coaxial illumination with 15V 150W lamp (bulb EFR 64634).
**should a lamp fail during an examination, the dual-port feature allows a quick change from one lamp to another resulting in an uninterrupted procedure.
- Maximum height from floor to objective lens – 46 inches.
- Minimum height – 34 inches.
- Green filter in selector knob.
- Small diameter base on roller castors with adjustable vertical post height.
- Spare lamp, dust cover and instruction manual.
- Video accessories available (see back page).

985



PRODUCT #	DESCRIPTION	PRICE
30985-OTS	Model 965 OTS Colposcope	\$6,995.00

Includes the following:

- Smooth five-step turret magnification changer
- 300mm working distance objective lens offering 3x, 4x, 7x, 11x, and 17x magnifications.
- Optional magnifications/objective lens combinations are available.
- Widefield 12.5x eyepieces with diopter locks and fold-down rubber eyecups for eyeglass wearers.
- Straight binocular head f=160.
- Fan cooled, **dual-port variable fiber optic coaxial illumination with 15V 150W lamp (bulb EFR 64634).
**should a lamp fail during an examination, the dual-port feature allows a quick change from one lamp to another resulting in an uninterrupted procedure.
- Maximum height from floor to objective lens – 54 inches. Minimum height – 32 inches
- Overall distance from the column to the objective lens – 33 inches.
- Green filter with selector knob.
- Spare lamp, dust cover and instruction manual.
- Video accessories available.
- Limited Lifetime Warranty

955



PRODUCT #	DESCRIPTION	PRICE
30955-5	Model 955 5-Step Colposcope	\$5,995.00

Includes the following:

- Smooth five-step turret magnification changer
- 300mm working distance objective lens offering 3x, 4x, 7x, 11x, and 17x magnifications.
- Optional magnifications/objective lens combinations are available.
- Widefield 12.5x eyepieces with diopter locks and fold-down rubber eyecups for eyeglass wearers.
- Straight binocular head f=160.
- Fan cooled, **dual-port variable fiber optic coaxial illumination with 15V 150W lamp (bulb EFR 64634).
**should a lamp fail during an examination, the dual-port feature allows a quick change from one lamp to another resulting in an uninterrupted procedure.
- Maximum height from floor to objective lens – 54 inches. Minimum height – 32 inches
- Overall distance from the column to the objective lens – 33 inches.
- Green filter with selector knob.
- Spare lamp, dust cover and instruction manual.
- Video accessories available.

SEILER 900 SERIES COLPOSCOPE ACCESSORIES

900 Series with Video Package		
935 Colposcope with video package	30935-3-V	\$8,975.00
955 Colposcope with video package	30955-5-V	\$10,175.00
985-OTS Colposcope with video package	30985-OTS-V	\$11,975.00

* Video Package Includes the Colposcope plus the following attachments:

Beam splitter 20/80	MC-A219S	\$1,495.00
Video camera adapter (1:1) f=65mm	MC-A218S	\$995.00
CCD color camera	SI3000021	\$995.00
19" Flat Panel color monitor	19MAGTV	\$695.00

900 Series with Digital Package

*935 Colposcope (Coolpix/Sony digital package)	30935-3-D	\$8,185.00
*955 Colposcope (Coolpix/Sony digital package)	30955-5-D	\$9,385.00
*985-OTS Colposcope (Coolpix/Sony digital pack)	30985-OTS-D	\$11,385.00

*Digital Package Includes the Colposcope plus the following attachments:

Beam splitter 50/50	MC-A219S	\$1,495.00
Digital camera adapter	MC-A292S	\$1,695.00

** Sony Camera is not included

Universal Camera Attachments

Universal digital camera adapter	SI81029990SX	\$2,495.00
Universal digital SLR camera adapter	SI8194101-SLR	\$1,795.00

**Please call us with your camera make and model before ordering a universal digital CA