

SEILER IQ®



Seiler IQ®	
Seiler IQ Floor Mount	IQ-100
Seiler IQ Wall Mount	IQ-200
Seiler IQ High Wall Mount	IQ-300
Seiler IQ Ceiling Mount	IQ-400
Seiler IQ Table Mount	IQ-500
(All of the above models are equipped with a fixed 45° head)	
With 0-180 Inclinable Head	IQ-180

Includes Balanced Optical Pod with:

- 0-180 degree inclined binocular head
- Counterbalanced pantographic arm
- DUAL LIGHT SOURCE = 50W Metal Halide with backup 15V 150W Halogen
- Automatic regulation of line voltage — 90V to 250V
- 10x Wide Field Eyepieces with Dioptic Locks
- Manual focusing control
- Infinity corrected/coaxial illumination
- Easily maneuverable optical head with tension control knobs and handles
- Objective lens 250mm
- 3 magnifications of 3.5x, 7x, and 11x with 250mm objective working distance lens

Video/Digital Accessories

- 202402-BS Beam Splitter
- 202402-VA Video Adapter with "C" Mount
- S1300002 CCD High Resolution Color Camera (with S video and composite output)
- 19MAGTV 19" Flat Panel LCD Color Monitor
- S18102990SX Universal camera adapter (For ALL Digital Point and Shoot Cameras)
- S18194101-SLR Universal camera adapter (For ALL DSLR Cameras)

Additional Accessories

- P6-200 Objective Lens 200mm
- P6-250 Objective Lens 250mm
- P6-300 Objective Lens 300mm
- P6-350 Objective Lens 350mm
- P6-400 Objective Lens 400mm
- 202-12.5X Eyepiece 12.5X (each)
- 202-16X Eyepiece 16X (each)
- 202402EXT6 600mm Extension Arm
- P6-FOC Fiber Optic Cable

