

SeilerScope



**SEILER**
PRECISION MICROSCOPES
A DIVISION OF SEILER INSTRUMENT COMPANY

SeilerScope

Achromat Binocular
Microscope For The
Physician's Office



Catalog #
SXS820

Description
SeilerScope Binocular Microscope

Includes

Achromat 4/0.10 (160)
Achromat 10/0.35 (160/0.17)
Achromat 40/0.65 (160/0.17R)
Achromat 100/1.25 (160/0.17R HI, oil)
L-shaped sturdy base with built-in 6V-20W transformer
Halogen illumination – 6V/20W Halogen bulb, #64250 (\$10.00)
Reversed four hole nosepiece
Coaxial coarse and fine focus mechanism
Focusing tension adjustment knob
Graduated, right hand double plated mechanical stage
50mm x 76mm - 0.1mm Vernier
2 - 10x Widefield Eyepieces (18)
Centerable Abbe condenser 1.25NA with iris
Rack and pinon mount
Inclined binocular tube ($f=1x$) (Seidentopf type)
Protective vinyl cover (D400-2301)



170 E. Kirkham Ave., St. Louis, MO 63119-1791
314-968-2282 • 800-489-2282 • Fax: 314-968-3601 • E-mail: micro@seilerinst.com