

Westlab II



SEILER
PRECISION MICROSCOPES
A DIVISION OF SEILER INSTRUMENT COMPANY

WESTLAB II

Binocular Microscope



CATALOG #	DESCRIPTION
SXSZ814	Westlab II Binocular Microscope

Includes:

41010/1	Main body, 6V 20W Tungston Halogen (easily switched from 110V to 220V). Reversed nosepiece. Sturdy 180 x 240mm base with rubber feet. Low position, coaxial coarse and calibrated fine focus control with graduations (2 microns per division). 40mm focusing range with tension adjustment and safety autofocus stop. Adjustable condenser. Built-in power supply with variable illumination control and separate on/off switch. Double plate specimen stage with right hand mechanical stage control. Includes dust cover(D400-2301), blue filter, immersion oil, spare bulb – 64250 ([®] 10. ⁰⁰) and instruction manual
410202	Seidentopf Binocular Head 30 degree inclined 360-degree rotation. Interpupillary distance 55 – 75mm
410301	Eyepiece WF 10x/18 (pair)
410401/S	Objective Achromat 4x/0.10
410402/S	Objective Achromat 10x/0.25
410404/S	Objective Achromat 40x/0.65, spring loaded
410406/S	Objective Achromat 100x/1.25, oil, spring loaded
410501	Abbe Brightfield condenser N/A 1.25, with integrated iris diaphragm and filter holder

ACCESSORIES:

410203	Trinocular head (Must subtract 410202 binocular head at \$313.00 when adding video)
410802	"C" Mount for Trinocular head
SI300002	CCD high resolution color camera, 460 lines
SI300035	14 inch color monitor
410901	Simple Pol Set with polarizer and analyzer
410306/S	Extra WF 10x/20 Eyepiece (price each)
410305	WF 10x/18.5 High Eyepoint Eyepiece (price each)



170 E. Kirkham Ave., St. Louis, MO 63119-1791
314-968-2282 • 800-489-2282 • Fax: 314-968-3601 • E-mail: micro@seilerinst.com