

MICROLUX III



 **SEILER**
PRECISION MICROSCOPES
A DIVISION OF SEILER INSTRUMENT COMPANY

Microlux III

Binocular Microscope



Catalog #

Description

SXSZ817

Microlux III Binocular Microscope

SXSZ817 DV

Microlux III with Dual View

Includes:

Semi-Plan achromat 4/ 0.10 160/-

Semi-Plan achromat 10/ 0.25 160/0.17R

Semi-Plan achromat 40/0.65 160/0.17R

Semi-Plan achromat Oil 100/1.25 160/0.17R

U-shaped wide stand with hand rest

45° inclined binocular tube, 360 degree rotatable with 387mm viewing height

Ball bearing mechanical stage, size 150mm x 135mm with movement of 76mm x 50mm

Abbe Condenser N.A. 1.25 with Aspheric lens iris diaphragm, with rack and pinion movement

Co-axial coarse and fine focusing system with auto focus stop

12V 20W Halogen, variable illumination with 90V-240V (bulb re-order number 64425)

UL Approved power supply

Blue filter

10x WF High eyepoint eyepieces 25mm (FN18)

Protective dust cover and instruction manual

Items subject to change.
Contact Seiler Instrument or check
our web site for more information.

06.05



www.seilerinst.com

SEILER

PRECISION MICROSCOPES

A DIVISION OF SEILER INSTRUMENT COMPANY

170 E. Kirkham Ave., St. Louis, MO 63119-1791

314-968-2282 • 800-489-2282 • fax: 314-968-3601 • email: micro@seilerinst.com