



SEILER
MEDICAL DIVISION
ENT MICROSCOPE



- 6 step turret magnification system (3 step turret available)
- LED Illumination + 150,000 LUX, 50,000 hour life
- Widefield eyepieces with dioptic locks:
10x, 12.5x 16x, and 20x
- True Coaxial Illumination="Shadow Free"
- Objective Lens: 175mm, 200mm, 250mm, 300mm,
and 400mm
OPTIONAL: Vario Focus Lens - 200mm-350mm
- Straight Head Standard (0-220° Inclinable Binocular
Head or Inclined Fixed 45° head available)
- Red Free Filter standard
- USB and 12V DC power support directly on optic pod
(charge Handycams, GoPro, DSLR, HD CCD)

ALPHA AIR SERIES

FOCUS ON PRECISION.



ACCESSORIES



ASSISTANT
TEACHING HEAD



0-220° BINOCULAR
HEAD



WIDE-FIELD
EYEPIECES



LASER FILTER



BEAMSPLITTER



EXTENDER
ROTATION DEVICE



IRIS DIAPHRAGM



GO PRO HERO



DSLR CAMERA
ADAPTER



HD LIVE VIDEO
CAMERA



S-VIDEO CAMERA



VARIO FOCUS

MOUNTING OPTIONS



FLOOR




HIGH WALL



WALL



CEILING



TABLETOP

